OP Technik *Surgical technique*

neon

superior spine system

Inhalt

Vorwort

Preface

4

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
5. Surgical technique – cervico-thoracic stabilization

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Preface

Preface

The spinal system neon is a modular system for posterior stabilization. Its application range extends from occipito-cervical via cervical to mid-thoracic applications.



neon implants are made from a titanium alloy in accordance with the standard DIN EN ISO 9001. Only the OC rods are made from pure titanium.

The present surgical technique describes the working steps necessary for the application of the neon system. The neon implants and the instruments are described in part 1 „The system“.

This surgical technique is not sufficient as the sole basis for the successful application of the spinal system neon. It is recommended to study and learn the operating technique with and from an experienced surgeon.


▲
Cervicale spondylotische
Myelopathie C3-C7
Cervical spondylotic
myelopathy C3-C7

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The small pictures at the bottom of the page show the chronological application of the instruments that are used as per the surgical steps on the double page. Pictures with instruments that had been used before are blue-coloured.

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The handling of the instrumentation and the introduction of the implants is shown by means of a vertebral column model. The general and repeat working steps are summarized in the following chapter.

- Assembly of the tissue protector and trocar
- Use of the cutter
- Depth determination
- Screw implantation
- Inserting the easy fits
- Inserting of neon spacers
- Use of connectors
- Removing the neon easy fits
- Fitting the cross link

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◀ Opening the cortex

Following identification of the entrance point, this is decortized with the neon awl (CS 1446) or a high speed cutter.

General assembly notes

1. Surgical technique – general assembly notes

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

◀ Assembly of the tissue protector and the trocar

The depth stop (CS 1440-2) with the locking nut (CS 1440-3) is set into the handle for tissue protector (CS 1440-1). The adapter (CS 1440-5, -6 or -8) or trocar (CS 1452-2 or CS 1452-4) is screwed from the front onto the tissue protector. The trocar is additionally secured by means of the locking nut CS 1452-1. The trocar sleeve (CS 1452-3 or CS 1452-5) is placed onto the trocar. Using the locking nut (CS 1440-3), the depth stop can be fixed in the required drilling depth.

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
◀ Using the cutter


The cutter (CS 1442-36) is used for the expansion of the cortex, if necessary.

CS 1440-1, -2, -3


CS 1440-5, -6, -8


CS 1452-2, -3, -4, -5

CS 1452-1

CS 1442-36

CS 1441-26, -35

General assembly notes

1. Surgical technique – general assembly notes

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
◀ Depth determination by means of the tissue protector

Drill under image intensifier control until the required depth is reached. The depth stop is set up to the AO-fitting of the drill or drill wire guide. The drilling depth can then be read off. The drill depth corresponds to the screw length.

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
◀ Depth determination by means of the measuring sleeve or the depth gauge


The measurement sleeve (CS 1440-10) is applied such that the guide wire or drill wire is situated within the groove of the measurement sleeve. The length of the screw can be read off from the measurement sleeve. Alternatively there is also the option to use the depth gauge (UT 1068-23) (without picture).

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

◀ Screw implantation

Without utilizing a guide wire the neon screws must be inserted by means of the screw driver, locking, cannulated (CS 1447). The screw is fitted by means of the knurled wheel of the instrument (see pictures).

Using guide wires either the screw driver, locking, cannulated (CS 1447) or the screw driver, rigid, cannulated (CS 1454) can be used. To detach the screw driver CS 1447 the guide wire must be removed before.

Utilizing a guide wire and a trocar sleeve as a tissue protector, the screw driver, rigid, cannulated (CS 1454) must be used.

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

◀ Inserting the easy fits

The neon easy fits, golden (CS 1411-1T), are inserted into the neon screws by means of the screw driver (CS 1461). The screw driver is self-holding so that the easy fit cannot become undone.

General assembly notes

1. Surgical technique – general assembly notes

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◀ Alignment of the easy fits

After screwing in the neon easy fits, these can be aligned, if necessary, such that the neon easy fits are positioned as vertically as possible relative to the frontal plane.

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◀ Use of neon spacers

In the event of a screw height offset, fitting a neon spacer will provide for height compensation. The neon spacer is fitted using the forceps for neon spacers (CS 1457). The spacer can also be retrofitted.

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◀ Fitting rod and connector

A bending of the rod in the frontal plane can largely be avoided by means of the connectors in 3 different lengths. In the sagittal plane, the rod should be bent in accordance with the physiological lordosis. Using the distraction forceps for neon connectors (CS 1459), the connector can also be pushed over an already contoured rod.

The connectors with the rod are placed on the easy fits. The grey end of the connectors must point to ventral.

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
◀ Inserting the inlay

The inlay (CS 1409-01T) must be placed in the connector by means of dissecting forceps. The inlay provides for an additional angle variability of 8° in all directions.

CS 1459

CS 1462-1

CS 1462-2

CS 1448

CS 1463-1

General assembly notes

1. Surgical technique – general assembly notes

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◀ Fixing the neon easy fits

In order to fit neon connector and screw, the nut for the neon easy fit (CS 1411-2T) is screwed on by means of the holder (CS 1462-1). The instrument is also self-holding so that the nut cannot drop out during insertion.

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◀ Tightening the easy fit hex nut

Using the box screw driver (CS 1462-2), the hex nut of the neon easy fit is firmly tightened. For holding the screw, the screw holder (CS 1448) can be used.

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◀ Removing the protruding part of the neon easy fits

The shear-off hex nut (CS 1463-1) is screwed onto the neon easy fit. For this purpose, the holder for the easy fit hex nut (CS 1462-1) can be used.

Important: Fit the hex nut with the thinner side facing downwards.

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▲ Removing the protruding part of the neon easy fits

The shear-off counter part (CS 1463-2) is fitted and the T-handled shear-off instrument (CS 1463-3) is slid into the counterpart. By fixing the counterpart and rotating the T-handle, the easy fit is sheared off.

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
◀ Tightening the neon easy fit hex nut

The shearing-off of the neon easy fit produces a smooth surface which is flush with the easy fit hex nut. Finally, retighten easy fit hex nut again by means of the box screw driver (CS 1462-2). The rod holder (CS 1458) is used as a counterpart. This will avoid forces acting onto the spinal column.

General assembly notes

1. Surgical technique – general assembly notes

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◀ Fitting the cross link

The hooks for the cross link (CS 1415-3T) are hooked into the rod. The hooks can be placed medially or laterally (picture) from the longitudinal support. After placing the X-plate for the cross link (CS 1415-1T/2T) the locking plate (CS 1415-5T) is attached by means of the easy fit hex nut using the box screw driver (CS 1462-2). The protruding threaded rods are broken off by bending with the box screw driver.

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The occipito-cervical (OC) stabilization can be carried out with two different implants, closed neon M plate or two neon OC rods. Two neon OC rods are easier to fit, in particular in anatomically difficult situations.

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2.1 Stabilization with neon M plate

The OC stabilization with the neon M plate (CS 1406T) is shown by the example of a C0-C2 instrumentation complete with transarticular screw union C1/C2 (see 3.1).

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Fitting the neon M plate for OC stabilization

At the beginning of the instrumentation, two neon screws are implanted transarticularly C1/C2 (see 3.1), and the easy fits are secured. The neon M plate rods are shortened to the correct length. The neon M plate is placed in position and pre-fixed to the neon screws by means of the connectors, using easy fit hex nuts.

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Drilling on the occiput and determining the screw length

The tissue protector (CS 1440-1/-2/-3/-5) is assembled and the drill depth adjusted. Using the drill Ø 2.6 mm (CS 1441-26), drilling is effected through the apertures of the neon M plate. The drilling depth indicated corresponds to the length of the thread. Therefore, the inserted screw must be longer by 2 mm (e.g. indicated drilling depth 8 mm, inserted screw 10 mm). Alternatively, the screw length can be determined by means of the depth gauge (UT 1068-23).


CS 1462-2

CS 1440-1, -2, -3, -5

CS 1441-26

UT 1068-23

CS 1970

Occipito-cervical stabilization

2. Surgical technique – occipito-cervical stabilization

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
◀ Use of neon toggles

neon toggles (CS 1417T) are used if bone thickness is less than 6 mm or if bone quality is poor. The neon M plate is removed again, and the bore is widened into a rectangular opening using a punch or bone cutter.

Important: The alignment of the punch mounting position must be in a transverse direction relative to the course of the plate. The neon toggles are inserted into the rectangular opening and turned through 90° for locking in that position.

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

◀ Fitting the neon M plate and fixing the toggles

The neon M plate is pushed over the toggles. At the same time the toggles must be held in their positions, e.g. by means of claws. The neon screws are used to fix the plate in position. First the locking plate, than the hex nut for the neon toggle is fitted.

Important: With its mandrels (see arrow in fig.), the locking plate must fit exactly into the plate aperture slots. This ensures rotational stability.

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◀ Fixing the neon toggle

The hex nut for the neon toggle is tightened by means of the flat spanner (CS 1970-X).

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◀ Implanting the cortical screws

If no neon toggle is required, cortical screws are used to fix the neon M plate. The cortical screws are implanted by means of the screw driver, hex 2.5 mm (CS 1328).

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◀ Removing the protruding ends

The protruding end of the neon toggle is removed by means of a wire cutter. To complete the assembly the protruding ends of the neon easy fits are removed.

◀ Komplette Montage der neon M-Platte


◀ Complete assembly of the neon M plate

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
2.2 Stabilization with neon OC rods

The use of neon OC rods is shown by the example of a C0-C3 instrumentation complete with transpedicular screw union C3 (see 4.2) and transarticular screw union C1/C2 (see 3.1).

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◀ Determining the form and length of the OC rod

Using the phantom rod (CS 1449), the required form and length of the neon OC rod are determined after the neon screws and easy fits have been set.

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
▲ Contouring the OC rod

Using the bender (CS 1460), the OC rod is contoured in the necessary planes and shortened to the correct length.

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◀ Prefixing the OC rod

The first fixation of the OC rod is effected by means of the connectors which are pushed onto the OC rod and attached to the neon screws using easy fits.

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



◀ Fixing the OC rod

Drilling in the occiput for cortical screws are carried out through the apertures of the OC rod by means of the drill 2.6 mm (CS 1441-26) and the tissue protector (CS 1440-1/-2/-3/-5). The cortical screws are implanted by means of the screw driver, hex 2.5 mm (CS 1328).

Alternatively, neon toggles (CS 1417T) can be used here also if bone thickness is less than 6 mm or if bone quality is poor (see 2.1).

Occipito-cervical stabilization

2. Surgical technique – occipito-cervical stabilization

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

◀ Assembly of the second OC rod

After assembling the second OC rod the instrumentation is finished by removing the protruding ends of the neon easy fits.

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3.1 Transarticular instrumentation and neon atlas claw

The following describes the transarticular C1/C2 instrumentation according to Magerl in connection with the neon atlas claw.

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

▲ Bore for transarticular screw

After opening the cortex the trocar tipped wire, Ø 1.5 mm (CS 1443-550) is introduced transarticularly. In order to keep the access small in the case of short length instrumentations, it is recommended to use the percutaneous instrumentation with the trocar system via two highly thoracic additional incisions. The trocar system comprises the tissue protector (CS 1440-1,-2,-3), trocar Ø 1.6 mm (CS 1452-2) and trocar sleeve (CS 1452-3), and it is inserted via a skin incision. The drill wire will then be inserted through the trocar system. At the depth stop the length of the screw can then be read off.

Atlantoaxial instrumentation

3. Surgical technique – atlantoaxial instrumentation

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▲ Use of the wire guide

If necessary, the wire guide (CS 1450) can be used if the rigid tissue protector does not allow for an optimum wire direction to be achieved.

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
◀ Wires introduced at both ends

Trocar sleeves situated in parallel before the screw insertion point.

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
◀ Implanting the neon screws

The self-drilling neon screw (CS 1400) of appropriate length is pushed over the trocar tipped wire (CS 1443-550), using the screw driver, rigid, cannulated (CS 1454), and implanted. Due to the drill tip, the screws should on no account be too long in order to avoid irritations of the pre-vertebral soft parts.

Important: When implanting the screw, hold the wire in its position.

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



◀ Fixing the neon easy fits

After implanting the transarticular neon screws, the wires are removed; the trocar sleeves remain in position until the operation has been completed. The neon easy fits are screwed into the screw heads.

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
◀ Fitting and inserting the neon atlas claw

The atlas hook for the neon atlas claw (CS 1416-1T) is shortened to the required length, the connector and the counter part (CS 1416-2T) are placed onto the end of the rod.

The atlas claw is hooked into the vertebral bow of the atlas and the connector complete with easy fit hex nut is prefixed on the screw. The atlas hook should not be fixed in a position further away than 16 mm lateral from the center line, in order to avoid injuries to the arteria vertebralis.

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◀ Positioning the counter part for the neon atlas claw

After tightening the easy fit hex nuts, the counter part is pushed under the caudal end of the atlas vertebral bow using the distraction forceps (CS 1465), and the neon atlas claw is closed.

The distraction forceps is supported by the connector and the counter part. The locking screw (arrow in fig.) of the counter part is tightened by means of the box screw driver (CS 1462-2). The protruding part of the neon easy fits is removed.

Atlantoaxial instrumentation

3. Surgical technique – atlantoaxial instrumentation

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
▲
Complete C1/C2 assembly
with neon atlas claw


▲
Instabile Denspseudarthrose
Unstable odontoid non-union

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

3.2 Isthmus and lateral mass

instrumentation

In cases where due to a variation in the course of the arteria vertebralis no safe transarticular screw placement C1/C2 without injury to the arteria vertebralis is possible, an alternative posterior C1/C2 instrumentation can be used. Self-tapping screws (CS 1401) are used which, in C1, are placed in the lateral mass and, in C2, in the isthmus.

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◀ Implanting of the neon screws in the isthmus

For isthmus screws, the trocar system with the short trocar is used (CS 1440-1, -2, -3, CS 1452-4, CS 1452-5).

After inserting the trocar system via a highly thoracic skin incision, the drill, 2.6 mm (CS 1441-26) is inserted.

The determination of the screw length is effected by means of the depth stop, the drill is removed and the guide wires (CS 1445-400) are inserted. After the trocar and tissue protector have been removed, the screws are inserted over the guide wire and through the trocar sleeves, using the screw driver, rigid, cannulated (CS 1454).

Atlantoaxial instrumentation

3. Surgical technique – atlantoaxial instrumentation

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◀ Implanting the neon screws into the lateral mass

For the bore in the lateral mass in C1, the tissue protector (CS 1440-1, -2, -3, -5) and the drill 2.6 mm (CS 1441-26) can be used. After removing the drill and introducing guide wires (CS 1445-400), the tissue protector is removed. The neon screws are implanted over the guide wire by means of the screw driver, rigid, cannulated (CS 1454).

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◀ Fixing the assembly

After screwing in the neon easy fits, the rods with the connectors are placed on laterally and fixed in position by means of the easy fit hex nuts.

CS 1440-1, -2, -3, -6

CS 1441-26

CS 1445-400

CS 1454


CS 1440-1, -2, -3, -5

CS 1441-26

Atlantoaxial instrumentation

3. Surgical technique – atlantoaxial instrumentation

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
▲ **Complete assembly with lateral mass and isthmus screws**

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3.3 Transpedicular and lateral mass instrumentation

A further alternative dorsal C1/C2 instrumentation is the implantation of self-tapping screws (CS 1401) in the lateral mass of C1 and, in C2, transpedicularly.

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◀ **Drilling for transpedicular screws**

The screws in C1 are placed in the lateral mass as described at 3.2. For the transpedicular screws, drilling is effected by means of the tissue protector (CS 1440-1, -2, -3, -5) and the drill Ø 2.6 mm (CS 1441-26). After removing the drill and introducing guide wires (CS 1445-400), the tissue protector is removed.

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◀ Implanting the neon screws

The transpedicular screws are applied via the guide wire, using the screw driver, locking, cannulated (CS 1447). Alternatively, the screw driver, rigid, cannulated (CS 1454) can be used.

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◀ Fixing the assembly

After screwing in the neon easy fits, the rods with the connectors are placed on medially and fixed in position by means of the easy fit hex nuts.

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

◀ Complete assembly with lateral mass and transpedicular screws.

Atlantoaxial instrumentation

3. Surgical technique – atlantoaxial instrumentation

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4.1 Lateral mass instrumentation

The subaxial instrumentation is shown by the example of a C3-C5 screw union of the lateral mass.

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◀ Drilling for lateral mass screws

For the bore in the lateral mass, the tissue protector (CS 1440-1, -2, -3, -5) and the drill Ø 2.6 mm (CS 1441-26) are used. After removing the drill and introducing guide wires (CS 1445-400), the tissue protector is removed.

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◀ Implanting the neon screws

The neon screws are implanted over the guide wire by means of the screw driver, locking, cannulated (CS 1447). As an alternative to the cannulated screws, non-cannulated screws Ø 3.5 mm (CS 1402) can also be inserted by means of the screw driver (CS 1447), but without any guide wires.

4

5

Subaxial instrumentations

4. Surgical technique – subaxial instrumentations

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◀ Determining the rod length

Following application of the neon easy fits, the contour of the required rod is determined using the phantom rod (CS 1449). The length of the rod is determined by means of the measuring caliper (CS 5788) (no picture).

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


◀ Fitting rod and connector

The neon rod complete with connectors is placed onto the screws via the neon easy fits, and fixed in position by means of the easy fit hex nuts. The placement of the rod complete with connectors is effected laterally or medially relative to the screw heads.

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
◀ Complete assembly C3 to C5 with lateral mass screws


CS 1449

CS 5788

Subaxial instrumentations

4. Surgical technique – subaxial instrumentations

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
◀ Drilling for transpedicular screws

For the transpedicular screws in the cervical spine, drilling is effected by means of the trocar system (CS 1440-1, -2, -3, CS 1452-4, CS 1452-5). The trocar system is introduced by lateral stab incision, and the drill Ø 2.6 mm (CS 1441-26) is inserted. After determining the screw length by means of the depth stop, the drill is removed and the guide wires (CS 1445-400) are inserted through the trocar. For the following steps, only the trocar sleeves remain in position.

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
◀ Implanting the neon screws


The screws are inserted over the guide wire and through the trocar sleeves, using the screw driver, rigid, cannulated (CS 1454). The implantation of the screws via the guide wire prevents a perforation of the pedicle or a screw breakout from the pedicle. This significantly reduces the risk of screw misplacements in comparison with non-cannulated systems.

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
◀ **Fixing the neon easy fits**

After implanting the transpedicular neon screws the guide wires and the trocar sleeves are removed. The neon easy fits are screwed into the screw heads. In the case of this instrumentation, the rod with the connectors is positioned medially. For completing the instrumentation, the easy fits are sheared off.

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◀ **Complete assembly with transpedicular screws C4 and C6**

5

Subaxial instrumentation

4. Surgical technique – subaxial instrumentations

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
Fitting in the cervico-thoracic area is demonstrated by means of transpedicular instrumentation from C6 to Th3.

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◀ Drilling in vertebra C6

For the transpedicular screws in C6, drilling is effected percutaneously by means of the trocar system (CS 1440-1, -2, -3, CS 1452-4, CS 1452-5). After inserting the drill Ø 2.6 mm (CS 1441-26), the length of the required rod is determined by means of the depth stop. The guide wire (CS 1445-400) is inserted into the pedicle through the trocar, and the tissue protector and the trocar are subsequently removed.

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


◀ Drilling in vertebra C7

For the bore in C7, the trocar system (CS 1440-1, -2, -3, CS 1452) is used. As drilling is not effected percutaneously, the trocar sleeve is not required. The guide wire (CS 1445-400) is inserted into the pedicle through the trocar, and the tissue protector and the trocar are subsequently removed. The view towards cranial shows the convergence angles that exist during drilling due to the pedicle course of C6 and C7.

Cervico-thoracic stabilization

5. Surgical technique – Cervico-thoracic stabilization

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◀ Implanting the neon screws

The implantation of the neon screws, Ø 4 mm (CS 1401) via the guide wire by means of the screw driver, locking, cannulated (CS 1447). Alternatively, the screw driver, rigid, cannulated (CS 1454) can be used.

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◀ Drilling in the thoracic spine

If screws with Ø 5 mm are used in the thoracic spine, the tissue protector (CS 1440-1, -2, -3, -8) and the drill Ø 3.5 mm (CS 1441-35) are used. When using screws with Ø 4 mm, the same instruments as on the cervical vertebral column are used. After determining the screw length by means of the depth stop, the drill is removed and the guide wires (CS 1445-400) are inserted through the tissue protector. Subsequently, the tissue protector is removed. The view towards cranial shows that the convergence angles, which arise during drilling due to the pedicle anatomy of Th2 and Th3, are significantly lesser than for C6/C7.

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




◀ Implanting the neon screws


In the thoracic vertebra, the neon screws, Ø 5 mm (CS 1403) are implanted via the horizontal guide wires, using the screw driver, rigid, cannulated (CS 1454).

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◀ Fitting

The neon easy fits are screwed into the screws. The rod with the connectors is positioned medially. For completing the instrumentation, the easy fits are sheared off.

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◀ Complete assembly in the cervico-thoracic transition

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
Implants	Product number
neon screw, self-drilling, cannulated, Ø 4 mm, length 34 mm	CS 1400-34T
neon screw, self-drilling, cannulated, Ø 4 mm, length 36 mm	CS 1400-36T
neon screw, self-drilling, cannulated, Ø 4 mm, length 38 mm	CS 1400-38T
neon screw, self-drilling, cannulated, Ø 4 mm, length 40 mm	CS 1400-40T
neon screw, self-drilling, cannulated, Ø 4 mm, length 42 mm	CS 1400-42T
neon screw, self-drilling, cannulated, Ø 4 mm, length 44 mm	CS 1400-44T
neon screw, self-drilling, cannulated, Ø 4 mm, length 46 mm	CS 1400-46T
neon screw, self-drilling, cannulated, Ø 4 mm, length 48 mm	CS 1400-48T
neon screw, self-drilling, cannulated, Ø 4 mm, length 50 mm	CS 1400-50T
neon screw, self-tapping, cannulated, Ø 4 mm, length 10 mm	CS 1401-10T
neon screw, self-tapping, cannulated, Ø 4 mm, length 12 mm	CS 1401-12T
neon screw, self-tapping, cannulated, Ø 4 mm, length 14 mm	CS 1401-14T
neon screw, self-tapping, cannulated, Ø 4 mm, length 16 mm	CS 1401-16T
neon screw, self-tapping, cannulated, Ø 4 mm, length 18 mm	CS 1401-18T
neon screw, self-tapping, cannulated, Ø 4 mm, length 20 mm	CS 1401-20T
neon screw, self-tapping, cannulated, Ø 4 mm, length 22 mm	CS 1401-22T
neon screw, self-tapping, cannulated, Ø 4 mm, length 24 mm	CS 1401-24T
neon screw, self-tapping, cannulated, Ø 4 mm, length 26 mm	CS 1401-26T
neon screw, self-tapping, cannulated, Ø 4 mm, length 28 mm	CS 1401-28T
neon screw, self-tapping, cannulated, Ø 4 mm, length 30 mm	CS 1401-30T
neon screw, self-tapping, cannulated, Ø 4 mm, length 32 mm	CS 1401-32T
neon screw, self-tapping, cannulated, Ø 4 mm, length 34 mm	CS 1401-34T
neon screw, self-tapping, cannulated, Ø 4 mm, length 36 mm	CS 1401-36T
neon screw, self-tapping, cannulated, Ø 4 mm, length 38 mm	CS 1401-38T
neon screw, self-tapping, Ø 3.5 mm, length 10 mm	CS 1402-10T
neon screw, self-tapping, Ø 3.5 mm, length 12 mm	CS 1402-12T
neon screw, self-tapping, Ø 3.5 mm, length 14 mm	CS 1402-14T
neon screw, self-tapping, Ø 3.5 mm, length 16 mm	CS 1402-16T
neon screw, self-tapping, Ø 3.5 mm, length 18 mm	CS 1402-18T
neon screw, self-tapping, Ø 3.5 mm, length 20 mm	CS 1402-20T
neon screw, self-tapping, cannulated, Ø 5 mm, length 30 mm	CS 1403-30T
neon screw, self-tapping, cannulated, Ø 5 mm, length 35 mm	CS 1403-35T
neon screw, self-tapping, cannulated, Ø 5 mm, length 40 mm	CS 1403-40T
neon screw, self-tapping, cannulated, Ø 5 mm, length 45 mm	CS 1403-45T
neon screw, self-tapping, cannulated, Ø 5 mm, length 50 mm	CS 1403-50T
Cortical screw, Ø 3.5 mm, length 8 mm	CS 1404-8T
Cortical screw, Ø 3.5 mm, length 10 mm	CS 1404-10T
Cortical screw, Ø 3.5 mm, length 12 mm	CS 1404-12T
Cortical screw, Ø 3.5 mm, length 14 mm	CS 1404-14T
Cortical screw, Ø 3.5 mm, length 16 mm	CS 1404-16T
Cortical screw, Ø 3.5 mm, length 18 mm	CS 1404-18T
neon M plate for OC stabilization	CS 1406T
neon OC rod	CS 1407T
neon rod, Ø 4.5 mm, length 20 mm	CS 1408-20T
neon rod, Ø 4.5 mm, length 25 mm	CS 1408-25T
neon rod, Ø 4.5 mm, length 30 mm	CS 1408-30T
neon rod, Ø 4.5 mm, length 35 mm	CS 1408-35T
neon rod, Ø 4.5 mm, length 40 mm	CS 1408-40T

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Implants	Product number
neon rod, Ø 4.5 mm, length 45 mm	CS 1408-45T
neon rod, Ø 4.5 mm, length 50 mm	CS 1408-50T
neon rod, Ø 4.5 mm, length 55 mm	CS 1408-55T
neon rod, Ø 4.5 mm, length 60 mm	CS 1408-60T
neon rod, Ø 4.5 mm, length 65 mm	CS 1408-65T
neon rod, Ø 4.5 mm, length 70 mm	CS 1408-70T
neon rod, Ø 4.5 mm, length 75 mm	CS 1408-75T
neon rod, Ø 4.5 mm, length 80 mm	CS 1408-80T
neon rod, Ø 4.5 mm, length 85 mm	CS 1408-85T
neon rod, Ø 4.5 mm, length 90 mm	CS 1408-90T
neon rod, Ø 4.5 mm, length 95 mm	CS 1408-95T
neon rod, Ø 4.5 mm, length 100 mm	CS 1408-100T
neon rod, Ø 4.5 mm, length 120 mm	CS 1408-120T
neon rod, Ø 4.5 mm, length 140 mm	CS 1408-140T
neon rod, Ø 4.5 mm, length 160 mm	CS 1408-160T
neon rod, Ø 4.5 mm, length 200 mm	CS 1408-200T
neon rod, Ø 4.5 mm, length 230 mm	CS 1408-230T
neon connector, polyaxial, green, length 15 mm	CS 1409-15T
neon connector, polyaxial, red, length 17 mm	CS 1409-17T
neon connector, polyaxial, blue, length 19 mm	CS 1409-19T
Inlay for neon connector, polyaxial	CS 1409-01T
neon locking screw for connector with length 19 mm	CS 1410T
neon easy fit, golden	CS 1411-3T
Hex nut for neon easy fit, golden	CS 1411-4T
neon spacer, red, height 2 mm	CS 1412-2T
neon spacer, green, height 4 mm	CS 1412-4T
neon cross link, length 50 mm, consisting of CS 1415-1T, -3T, -5T	CS 1415-0-50T
neon cross link, length 65 mm, consisting of CS 1415-2T, -3T, -5T	CS 1415-0-65T
X-plate for neon cross link, length 50 mm	CS 1415-1T
X-plate for neon cross link, length 65 mm	CS 1415-2T
Hook for neon cross link	CS 1415-3T
Locking plate for neon cross link	CS 1415-5T
neon atlas claw consisting of CS 1416-1T, -2T, -3T	CS 1416T
Atlas hook for neon atlas claw	CS 1416-1T
Counter part for neon atlas claw	CS 1416-2T
Locking screw for neon atlas claw	CS 1416-3T
neon toggle consisting of CS 1417-1T, -2T, -3T	CS 1417T
neon toggle	CS 1417-1T
Locking plate for neon toggle	CS 1417-2T
Hex nut for neon toggle	CS 1417-3T

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Instruments	Product number
Screw driver, hex 2.5 mm	CS 1328
Handle for tissue protector	CS 1440-1
Depth stop for tissue protector	CS 1440-2
Locking nut for tissue protector	CS 1440-3
Adapter for tissue protector, Ø 2.6 mm	CS 1440-5
Adapter for tissue protector, Ø 2.6 mm, blunt	CS 1440-6
Adapter for tissue protector, Ø 3.5 mm	CS 1440-8
Measurement sleeve for determining the depth	CS 1440-10
AO drill bit, Ø 2.6 mm	CS 1441-26
AO drill bit, Ø 3.5 mm	CS 1441-35
Cutter, Ø 3.6 mm, cannulated	CS 1442-36
Trocar tipped wire, Ø 1.5 mm, length 550 mm	CS 1443-550
Guide wire, blunt, Ø 1.5 mm, length 400 mm	CS 1445-400
neon awl	CS 1446
Screw driver, locking, cannulated, const. of CS 1447-1, CS 1447-2	CS 1447
Handle for screw driver, locking	CS 1447-1
Attachment for screw driver, locking	CS 1447-2
Screw holder	CS 1448
neon phantom rod, length 250 mm	CS 1449
Wire guide	CS 1450
Locking nut for trocars	CS 1452-1
Trocar, cannulated, Ø 1.6 mm, length 173 mm	CS 1452-2
Trocar sleeve for CS 1452-2	CS 1452-3
Trocar, cannulated, Ø 2.7 mm, length 137 mm	CS 1452-4
Trocar sleeve for CS 1452-4	CS 1452-5
Screw driver, rigid, cannulated	CS 1454
Forceps for neon spacer	CS 1457
Rod holder, Ø 4.5 mm	CS 1458
Distraction forceps for neon connector	CS 1459
Bender for neon rod and neon OC implants	CS 1460
Screw driver for neon easy fit	CS 1461
Holder for neon easy fit hex nut	CS 1462-1
Box screw driver for neon easy fit hex nut	CS 1462-2
Shear-off hex nut	CS 1463-1
Shear-off counterpart	CS 1463-2
T-handled shear-off instrument	CS 1463-3
Compression forceps	CS 1464
Distraction forceps	CS 1465
Flat spanner, length 120 mm	CS 1970-00-12
Flat spanner, angle 45°, length 120 mm	CS 1970-45-12
Measuring caliper	CS 5788
Depth gauge, length 230 mm	UT 1068-23

ulrich GmbH & Co. KG

Buchbrunnenweg 12

89081 Ulm

Germany

Telefon/Phone +49 (0)731 9654-110/225

Telefax/Fax +49 (0)731 9654-2702

e-mail spine@ulrich-ulm.de

internet www.ulrich-ulm.de